

- The first rechargeable Coselgi BTE hearing aid
- Up to 24 hours of usage inclusive eight hours of streaming per charge
- 4 hours for full charge, 30 minutes fast charge for 4 hours of usage
- Easy adjustment to individual requirements

## REAL-LIFE HEARING


Since 1961, Coselgi has been delivering the advanced hearing technology people need to accurately reflect the world around them. We understand that real life doesn't happen in a lab or office, and we work hard to combine our technologies with professional expertise to deliver real-life sound for real-life hearing.

[www.coselgi.com](http://www.coselgi.com)



COSELGI is a trademark of WS Audiology A/S

V1 / 10.2022



## REAL-LIFE SOUND FOR REAL-LIFE HEARING



Mojo hearing aids use advanced technology to produce real-life sound, accurately representing the world around you.



## MOJO HEARING AIDS

Hearing is highly personal, and when you want to hear what's important to you, having an accurate representation of your world is essential. Mojo hearing aids with bluetooth technology use advanced technology that combines data from how you hear and how you use your hearing aids to deliver sound that's tailored to you. With a broad range of models, including rechargeable options, Mojo hearing aids are the best way to get real-life sound for real-life hearing.



# REAL-LIFE SOUND, TAILORED TO YOU

The advanced technology in Mojo hearing aids works to produce real-life sound right from the first fit. As you continue to use them day-to-day, intelligent automation analyzes the sounds around you and continuously selects the signal processing most appropriate for the situation you are in.

When you use the sound optimization tool, the Real-Life Preference Log remembers your personal preferences and personal programs you have stored. When shared with your hearing care professional, this enables them to fine-tune your hearing aids based on how you actually use them – resulting in a hearing experience best suited for your world around you.



# BTE R D – OUR NEW RECHARGEABLE MODEL

While developing this hearing aid we were focused on real-life sound for real-life hearing and gave special interest to operating convenience. No more fumbling to change the battery and easy control of the hearing aid with the preference control and separat program button. You can get an indication of the battery status simply by looking at the integrated LED in the hearing aid.



BTE R D with  
charger

# A MODEL TO SUIT EVERY NEED

Coselgi Mojo hearing aids come in six behind-the-ear models and three in-the-ear models. Your audiologist will recommend which model best meets your needs. When it comes to color though, you call the shots!



- Autumn Beige
- Silver Grey
- Titanium Grey
- Tech Black
- Golden Brown
- Honey Blonde
- Chestnut Brown



Relative size  
to iPhone 13

Sounds from smartphones, smart speakers, TVs and tablets are a normal part of life these days. Using the Coselgi Mojo App, you can take full advantage of this by seamlessly connecting to iOS and Android (Android 10 and higher) devices to stream high quality audio straight to your hearing aids - and make individual adjustments to your sound preferences. The Coselgi Mojo app is compatible with Mojo hearing aids with bluetooth technology.



Download the  
Coselgi Apps now!